

CeeMax™ 802ss Table Top Slit Lamp

Slit Illumination

Type	Galilean-Type
Magnification Ratio	6X, 10X, 16X, 25X, 40X
Objective	F = 108mm
Eyepieces	12.5X
Angle between eyepieces	13°Converging optics
Field of view	35mm, 23mm, 14mm, 9mm, 6mm
Pupillary Adjustment	51mm~78.5mm
Diopter	-6 to+6 Diopter
Slit Width(mm)	Continuously variable from 0 to 13(at 13mm, slit becomes a circle)
Slit Length(mm)	Continuously variable from 1 to 13
Aperture Diameters(mm)	13, 10, 5, 3, 2, 1, 0.2
Slit Angle	0°~180°Continuously variable
Slit Inclination	4 steps: 5°, 10°, 15°, 20°
Filters	Heat absorbing, UV, Red-free, Cobalt Blue
Lamp	10W LED
Illumination Type	From top(Haag-Streit type)

Chin Rest

Vertical movement	70mm
Fixation Target	Green LED

Base

Longitudinal movement	90mm
Lateral movement	100mm
Fine Base movement	15mm
Vertical movement	30m

Demension & Weight

Package Dimension	760mmx450mmx370mm(2)
Total Weight	24Kg(53 lb)

Power Supply

Input
Frequency

110V/230V AC
60Hz/50Hz



- Ready for digital photography and video
- Ready for Goldmann tonometer(Haag-Streit 900 type)
- One year factory warranty(contact CeeMax Ophthalmic)

Contact: **CeeMax™ Ophthalmic**

17100 Gillette Ave. Unit 113, Irvine, CA 92614 USA

Tel: 949 885 8875

email: info@ceemaxusa.com